

Order form Sigvaris Roundknit Custom Made Stockings

Leg measurements

SIGVARIS
GROUP

Country / Distributor

Contact person

Patient name

Sex ☐ female ☐ male

Order No.

Your Reference No.

Re-order No.

Order date

Remarks

SIGVARIS AG

9014 St. Gallen

Switzerland

Tel. +41 71 272 40 00

Fax +41 71 272 40 01

custom-made@sigvaris.com

www.sigvaris.com

sigvaris

Product code	New product name	Old product name
SSE	Style SEMITRSPARENT	MAGIC
ECF	Essential COMFORTABLE	COMFORT
ETH	Essential THERMOREGULATING	COTTON
ECL	Essential CLASSICAL	TOP FINE
TRA	Specialities TRADITIONAL	TRADITIONAL

Quality:

CCL:

CCL:

SSE

1 2

☐ ☐

ECF

1 2

☐ ☐

ETH

1 2 3

☐ ☐ ☐

ECL

1 2

☐ ☐

TRA

2 3 4

☐ ☐ ☐

☐ For Specialities TRADITIONAL compression up to the cG level (without soft-part)

Models:

A-D Calf

A-D with grip top

A-F Half-thigh

A-F Half-thigh with grip top*

A-G Thigh

A-G with grip top

A-G Thigh with waist attachment:

☐ left ☐ right

A-T Pantyhose

A-T Bodyform

A-T Maternity

A-T for men

A-T One leg Pantyhose

Colors:

beige

caramel

crispa

nature

savanna

black

skin

grey

marine

TREND COLORS: salted caramel

blueberry dream

fresh basil

hot chili

Pieces

Pairs

Grip tops:

☐ Sensinova grip top neutral design 5.0 cm

☐ Sensinova grip top feminine design 5.5 cm (only available for SSE)

☐ Knobbed grip top 5.0 cm

☐ Narrow knobbed grip top 2.5 cm (for A-D calf stockings available only)

☐ Micro knobbed grip top 5.0 cm (only available for ETH A-G thigh stockings in CCL1 and 2)

* for SSE only available in CCL 2

1 / 2

Order form Sigvaris Roundknit Custom Made Stockings

Leg measurements

Taking measurements

Please use a separate order form for each model.

When to measure

The measurements should be taken in the morning, on the edema-free leg, in a standing position.
The measurements should be made directly on the skin. Both legs must be measured, if garments for both legs are required.


How should be measured

Circumference: The tape measure should be tight but not cut in when measuring.
Length: Measure the length (with the exception of ℓT) on the inside of the extended leg in a perpendicular position.

NOTE

We cannot be held responsible for incorrect measurements given. Custom-made compression stockings cannot be returned.
Product liability cannot be accepted if the label has been removed.

Leg measurements:

Length Left leg	Circumference Left leg		Circumference Right leg	Length Right leg
ℓT	<i>Length of panty-part KT:</i>		cT	ℓT
ℓH	<i>KT length in front:</i>		cH	ℓH
	<i>KT length at the back:</i>			
ℓG	cG		cG	ℓG
ℓF	cF		cF	ℓF
ℓE	cE		cE	ℓE
ℓD	cD		cD	ℓD
ℓC	cC		cC	ℓC
ℓB^1	cB^1		cB^1	ℓB^1
ℓB	cB		cB	ℓB
	cY	cY		
	cA	cA		
	ℓA Foot length without toe	ℓA Foot length without toe		
	ℓZ Foot length incl. toe	ℓZ Foot length incl. toe		

Toe:

☐ Open toe

☐ Closed toe